MAILLI

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朔盈
Extreme

EN

For an electronic copy of the MAILI Instructions for Use please refer to: www.sinclairpharma.com/eifu and input the KEYCODE: MLX.COR. If you require a paper copy then please contact: quality@sinclairpharma.com or +44 (0) 1244 625150.

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Sterilised using moist heat (syringes)
经蒸汽灭菌 (注射器)

Sterilised using irradiation (needles)
经辐射灭菌 (注射针)

European Community Representative
欧盟授权代表

Sterile fluid path - Moist heat (the content of the syringe (gel) is steril)
包含经蒸汽灭菌的无菌液路 (注射器的内容物 (凝胶) 无菌)

CE-mark (needles)
CE 标志 (注射针)

CE-mark (syringes)
CE 标志 (注射器)

Description

MAILI EXTREME is a sterile, transparent and resorbing gel of cross-linked hyaluronic acid of biofermentative origin. The product is supplied in plastic syringe and it is intended for single use only.

The volume of each syringe is 1 mL for **MAILI EXTREME**. Each box contains four graduated pre-filled syringes of **MAILI EXTREME**, an instruction leaflet, a set of traceability labels and four single use sterile needles of 27G 1/2". For each syringe, one label is to be given to the patient, the other one is to be added to the patient record.

The content of **MAILI EXTREME** syringes is sterilised by moist heat. The 27G 1/2" needles are sterilised by radiation. The intended purpose of **MAILI EXTREME** is to modify the anatomy and/or to alleviate for a physiological process at the level of the face.

Composition

MAILI EXTREME

- Cross-linked* hyaluronic acid: 24 mg/mL
- Lidocaine hydrochloride monohydrate: 3 mg/mL
- Phosphate buffer: q.s. 1 mL
- One syringe contains 1 mL of **MAILI EXTREME**

This product consists of a pre-filled syringe, a disposable sterile stainless steel injection needle and a gel encapsulated in the syringe. The gel is composed of crosslinked sodium hyaluronate, lidocaine hydrochloride and phosphate buffer, and sodium hyaluronate is prepared by microbial fermentation. The labeled concentration of crosslinked sodium hyaluronate is 24mg/ml and the labeled concentration of lidocaine hydrochloride is 3 mg/ml. The syringe pre-filled with gel is sterilized by moist heat, and the injection needle is sterilized by gamma radiation. This product is for single use. The shelf life of the product is 2 years.

*The hyaluronic acid has been cross-linked with BDDE.

Applicable population

MaiLi Extreme is intended for adults 18 years and older.

Intended uses

MaiLi Extreme is indicated to inject into facial suprapreosteal zone, fill the jaw area to improve the mandibular contour of patients with mild to moderate mandibular retraction.

Mode of action

MAILI EXTREME creates a volume in the tissue that corrects facial anatomy.

MAILI EXTREME is intended to be wholly resorbed over time, which progressively reduces the volumizer effect of the product at the injection site. The kinetic of this absorption depends on several factors among which the injected quantity, the depth of injection and the patient metabolism.

In a clinical trial in China, a total of 159 subjects received submaxillary injection of MaiLi Extreme, most of whom were assessed as improving at the primary endpoint (6 months after the last treatment), with the majority of subjects continuing to have an effect at 12 months after the last treatment.

Contra-indications

- Do not inject patients having known risk of hypersensitivity to hyaluronic acid or lidocaine or amide local anaesthetics.
- Do not inject patient having history of severe allergy or anaphylactic shock.
- Do not inject patients suffering from autoimmune disease.
- Do not inject patients suffering from porphyria.
- Do not inject patients having cutaneous disorders, inflammation or infection (herpes, acne, etc.) at the treatment site or nearby.
- Do not inject patients for whom the medical history shows a sensitivity that could lead to a reaction to the treatment.

- Do not use in patients with bleeding disorders or in patients who are undergoing treatment with thrombolytics or anticoagulants.
- Do not inject patients with a tendency to form keloids, hypertrophic scars or any other healing disorders.
- Do not use during pregnancy or breast-feeding.
- Do not use for children.
- Do not inject patients with mandibular deformity and clinical symptoms

Warnings

- Do not inject into blood vessels (intra-vascular) or directly near blood vessels (for avoiding vessels compression). An unintended intravascular administration can cause high blood concentration and acute central nervous system and cardiovascular toxic symptoms.
- Do not inject in area other than the face (e.g. hands, body).
- Do not inject **MAILI EXTREME** in the periorbital area (eyelids, palpebromalar groove, crow's feet) and glabellar region.
- The injection of **MAILI EXTREME** into the subcutaneous fat tissue or into the suprapreosteal zone is reserved to specialists specifically trained in these injection techniques.
- Do not inject intramuscularly. There is no available clinical data about injection of **MAILI EXTREME** into an area which has already been treated with another brand of resorbing product. It is therefore recommended not to inject in areas in which another brand of resorbing product is not still fully resorbed.
- Do not inject jointly with another brand of resorbing product in the same area.
- Do not inject jointly with a permanent product.
- Do not mix with other products.
- Do not re-sterilise the product (the syringe and/or the needle).
- Do not inject more than 4mL in the first injection. If the best aesthetic results are not achieved, a supplementary injection can be given 4 weeks after the first injection, and the supplementary injection should not exceed 2ml. Do not inject more than 20 mL per year per patient.
- Do not overcorrect.
- Do not use after expiry date or if the packaging is damaged or the syringe opened.
- Do not re-use (risk of contamination). The product is intended for single use only.
- Do not use a syringe for several patients (risk of cross-contamination).
- Do not use simultaneously with chemical peels, dermabrasion, laser treatment or radiofrequency treatment.

Precautions for use

- Patients must be questioned regarding their medical history and must be informed about expected outcome, warnings, contra-indications, precautions and potential adverse events.
- The quantity to inject depends on the area to treat and the expectation of the patients. Patients with unrealistic expectation must not be treated.
- Aseptic techniques and standards practices are to be observed.
- The content of the syringe is sterile until opening. External surface of the syringe is not sterile and should not be used in a controlled environment such as operating room. If needed, the surface of the syringe must be cleaned up in respect of the aseptic techniques in place.
- The patients on medication that affects platelet functions (aspirin or non-inflammatory drugs) must be warned of the potential increased risks of haematomas and bleeding during injection.
- Special caution must be exercised when treating areas in close proximity to vulnerable structures

such as vessels and nerves.

- Advise the patients to avoid any extended exposure to the sun, UV rays and temperatures below 0°C, as well as sauna or hammam sessions during the two weeks following the injection procedure.
- Recommend the patients not use cosmetic products (e.g. makeup, skincare products) during the 12 hours following the injection treatment.
- Recommend the patients to avoid applying intense pressure or massaging the treatment site for a few days following the injection.
- The healthcare professional must consider the fact that the product contains lidocaine. For normal healthy adults, the maximum total dose of lidocaine should not exceed 200 mg per session. When using concurrently (case of topical administration), the total administered dose of lidocaine should be considered. The concomitant use of other local anaesthetic agents should also be considered since the systemic toxic effects may be additive. Care should be taken for patients with congenital methemoglobinemia and patients who are receiving concomitant treatment with methaemoglobin-inducing agents. In the case of athletes, the presence of lidocaine may produce a positive result in anti-doping tests.

Assembly of needle to syringe

- Use of the supplied 27G 1/2" sterile needle is recommended. Use of another needle could not ensure the correct fitting to the syringe and/or the correct delivery of the product.
- Remove the tip cap of the syringe.
- Grasp the hub.
- Press and turn the needle clockwise on the syringe until it is tight.
- Make sure that the hub of the needle is as tightly fixed to the bottom of the luer-lock as possible and remove the shield of the needle just before injection.
- Failure to comply with these precautions could cause disengagement of the needle and/or product leakage at luer-lock level.

Interactions

- Treatment with MAILI EXTREME in combination with drugs and other devices has not been clinically evaluated.
- Pre-treatment with a lidocaine solution is not required since the product already contains lidocaine.

Procedure of administration

Treatment with MAILI EXTREME is only intended to be administrated by an authorized healthcare professional in accordance with the local legislation.

The healthcare professional shall have a deep knowledge of the anatomy of the area to be treated and shall be used to perform similar intervention (i.e. cross-linked hyaluronic acid injection in the face area). MAILI EXTREME is to be used as supplied. Modification or use of the product outside of the instructions for use may adversely impact the sterility and/or the performance of the device.

A correct injection technique is important for the clinical outcome of the treatment. The injection technique varies depending on the treatment site, the depth of injection and the gel quantity to be injected.

- 1) Before starting treatment, the patient should be informed of the product's indications, contra-indications, precautions, warnings, expected results and potential adverse events and questioned on her/his medical history.
- 2) The injected area must be disinfected with an appropriate antiseptic solution.
- 3) Assemble the supplied sterile needle and syringe.
- 4) Do not apply extensive pressure to the syringe at any time. If such resistance is encountered, remove

and replace the needle and check for function.

5) MAILI EXTREME be administrated slowly in the tissues. The injection volume varies depending on the correction required. Do not overcorrect.

MAILI EXTREME is indicated for injection into the subcutaneous fat tissue or into the supraperiosteic zone.

6) For optimal clinical results, gentle massage by the healthcare professional of the treated area is recommended to ensure uniform distribution of the gel.

Adverse events

Healthcare professional shall inform the patients about potential adverse events related to the injection procedure and/or the use of the product MAILI EXTREME, occurring immediately after injection or after a delay.

Based on literature review on injectable hyaluronic acid gels and in agreement with the clinical data obtained with MAILI devices (clinical study 17E1176), based on the clinical trial data in domestic and globally, these adverse events include but are not limited to:

- Events which are naturally resolve within one week in most cases:
- Injection-related events and/or inflammatory reactions such as bleeding, ecchymosis, erythema, haematoma, skin redness, bruising, swelling, oedema and infection which may be associated with local pain or itching, occurring after injection.
- Sensitivity at the injection site.
- Hardness, lump or nodule at the injection site.
- Skin coloration or discoloration at the injection site.

- Events which are delayed after the injection:
- Immediate or delayed hypersensitivity to hyaluronic acid and/or to lidocaine.
- Nodule, abscesses or granuloma.
- Vascular compromise may occur due to inadvertent intravascular injection or as a result of vascular compression associated with implantation of any injectable soft tissue filler. This may manifest as blanching, discoloration, necrosis or ulceration at the implant site or in the area supplied by the blood vessels affected; or rarely as ischemic events in other organs due to embolization. Isolated rare cases of ischemic events affecting the eye leading to visual loss, and the brain resulting in cerebral infarction, following facial aesthetic treatments have been reported. These rare cases of vascular events are mostly reported in glabella, in and around the nose, in forehead and in periorbital area.
- Infection or reactivation of a previous infection.
- Displacement of the gel.

Inflammatory reactions which persist for more than one week, or any other adverse event which develops, must be reported by the patients to their healthcare professionals. In this case, if required, the healthcare professionals should use an appropriate treatment. Any adverse events and complications must be reported to the local representative of KYLANE LABORATOIRES.

Waste disposal

At the end of the session, opened and used syringe as well as all injection materials (e.g. needle) shall be discarded in respect of the medical practices and in accordance with the applicable regulations on waste disposal (e.g. specific waste container).

Storage conditions

- Shelf life of the product is indicated on each syringe and product must not be used after this date.
- Store between 2°C and 25°C away from direct sunlight. Keep container in the outer carton.
- Make sure there is no visible sign of damage to the packaging before use.

Manufactured by

Kylane Laboratoires SA
Chemin Pré-Fleuri 1-3
CH-1228 Plan-Les-Ouates
Switzerland

Product Name: Crosslinked Hyaluronan(HA) injectable gel with Lidocaine
Product Model: Maili Extreme
Specifications: 1ml/Syringe, 4 Syringes /Box
Medical Device Registration Certificate Number: 国械注进20253130026
Product Technical Requirements Number: 国械注进20253130026
Manufacturing Date: See original label
Use by date: See original label
Registrant name: Kylane Laboratoires SA
Address of Registrant: Chemin Pré-Fleuri 1-3, CH-1228 Plan-les-Ouates, Switzerland
Registrant Contact Information: +41 22 588 96 11
Manufacturing Site Name: SOTHEMA
Manufacturing Site / Plant Address: BP. N°1 Industrial area Bouskoura, 27182 Casablanca, Morocco
After-sales service unit /Agent Name: Sinclair (Hangzhou) Supply Chain Management Co., Ltd
After-sales service unit /Agent Address: Room 213, 2nd Floor, Building 7, No.161, 17th Street, Baiyang Street, Qiantang District, Hangzhou City, Zhejiang Province.
After-sales service unit /Agent Contact Information:
E-mail: Customerservice@SinclairAesthetics.com
Phone: +4008206181
Country of Origin: Morocco
The Revision date of IFU: 2025-03

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Sinclair France SAS